

**Model number:** PTCE-III Patient Goggle



- Stainless steel
- Protection for **all** laser and IPL wavelengths
- Maneuverable nose bridge flips 180 degrees
- Adjustable pupillary distance
- Detachable eyecups - can be used individually
- Silicone-wrapped eyecups
- Can be cold sterilized or autoclaved
- Silicone strap and hard case included
- Custom-engraved logo on case available

**Protection Ratings:**

315-1400nm (OD 7+)

1400-11000nm (OD 5+)

**CE:**

315-1400 D LB6 + IR LB8 + M LB8Y

1400-11000 D LB4 + I LB5 INVO CE



### Eye Shield Specifications

